



ADENSYS
PRODUCT CATALOGUE.
2022/2023

Designed to perform better Home of dental components



The ADENSYS Story

Founded in 2012, ADENSYS is an engineering and manufacturing firm based in the small Alpen town of Triesen, in Liechtenstein. The company was founded by Frits Diepen and Helmut Gruber, both highly qualified development engineers and managers with many years of previous experience with medium-to-large dental firms situated in the surrounding area. ADENSYS has a branch office (ADENSYS GmbH) located in Buxheim, in southern Germany, and it is here that most assembly and production take place. Being situated between Austria, Switzerland, Germany and Italy, ADENSYS is in a unique position to take advantage of the substantial local competences in micro-machining and dental technology.

Our Passion

We have passion for all things dental. We have a huge amount of knowledge and experience related to the dental industry and we feel that our know-how can help you better understand your products and your customer. We want to work with you, our customers and suppliers, to help develop your business. We believe that we can succeed only if you succeed.



Helmut GruberADENSYS founder

Frits DiepenADENSYS founder

Version: 220512 v4

Section 1 Parts for contra-angle/ attachment handpieces

Section 2 Parts for high-speed turbine/handpieces

Section 3Complete product overview

Section 4Additional information

About this book.

This catalogue is designed to present all of the products that are currently available from ADENSYS.

Our products are aimed at professionals that are involved in the repair of dental handpieces and, given this, we understand that different people need different points of reference when searching for parts. For this reason, we have presented our products in two ways:

i) by handpiece model and

ii) as collections by category.

We hope this will make it convenient to find the needed parts. Below, you will find some information about the different sections of this catalogue.

On the cover page of each section, you will find a table of contents to help you navigate through the pages.

Please note that you can order all products in this catalog at **www.adensys.net**.


All default prices at our website are street (i.e. dentist) prices. If you are a dental handpiece repair professional, please contact us by e-mail at **info@adensys.net** to see if you qualify for a wholesale discount on our listed prices.

About this book.

Model Pages

Sections 1 and 2 of this catalogue show details of the parts available per individual handpiece model. Please see the explanation below to understand how the information on these pages should be interpreted or understood.

If anything is unclear here, we encourage you to contact us at info@adensys.net so that we can evaluate any possible changes or improvements.




About this book.

Bearing pages

All individual ADENSYS parts are listed in section 3 of this catalogue. As we describe in Section 4, ball bearings are complex products themselves and can have many mechanical properties that need to be specified. Please see the guide below as to how ADENSYS bearings are described.

In an effort be transparent about our suppliers for our available bearing sizes, we have indicated the exact manufacturer for our Made-in-Germany bearings. Given our location in Germany and our positioning as a high-quality manufacturer, we try to use the highest-quality premium bearings in all of our applications. We do realize, however, that economics do play a critical role in many repair decisions, so we also offer a more competitively-priced alternative for some of our bearing models. In the explanation chart below, you will see that we have indicated the bearing manufacturer for the German-sourced products.



About this book.

Glossary

We use standard nomenclature to describe many of the product types and mechancial properties described in this book, and we know that there are often different expressions used for these things around the world. Forthis reason, please see the following definitions:

High-speed handpieces All air-driven or turbine handpieces, which turn a turbine-style rotor to

deliver torque to the cutting element/bur

Attachment handpieces All motor-driven handpieces, which feature an e-type attachment or

coupling (see more information on the e-type attachemnt below). These handpieces are often known as "contraangles" because of the bend in the middle of the handpiece.-style rotor to deliver torque to

the cutting element/bur

Head gear rotor Found in the head housing of attachment handpieces, this is the rotor

assembly that holds the bur or cutting instrument. This assembly is

sometimes also referred to as "cartridge".

Intermediate shaft Found in the neck of attachment handpieces, this rotor typically

consists of a long shaft and will mesh with the head gear rotor.

Complete Turbine Found in the head housing of high-speed handpieces, this is the rotor

assembly that holds the bur or cutting instrument.

Combo A sub-assembly, which consists of a spindle and impeller. This product

is dynamically balanced prior to delivery. Fitting bearings can be

assembled on this product to create a complete turbine.

Premium turbinesWithin our selection of complete turbines, these products are

assembled with top-of-the-Line bearings, equivalent to those use by

the original equipment manufacturer (OEM).

Standard turbines Within our selection of complete turbines, these products are

assembled with superior-quality ceramic bearings, usually from

Japanese maufacturers

Ceramic bearing

This term refers to the fact that we use bearings with ceramic balls. The

outer ring and inner ring are made from stainless steel. You will find

more information on bearings in section 4 of this book.

e-type connection Described the standard connection between typical electric or air

motors for dental application. The connection properties are described

in the norm ISO 3964

REF numberWherever "REF" is indicated in this book, it denotes the article number

from the original manufacturer of the handpiece or medical device.

Please note that ADENSYS products products will never be a 100%

copy of the original product, but it is 100% compatible


Parts for Attachment Handpieces.

On the following pages you will find for different manufacturers and their contra-angle handpieces all suitable spare parts for a complete repair.


KaVo MASTERmatic Lux M25L	10
KaVo MASTERmatic Lux M05L	11
KaVo GENTLEpower 25LP/25LPA/25LPR	12
KaVo EXPERTmatic E25L/E25C	13
KaVo INTRAmatic 25LH/25LHA/25LHN/25 LHC	14
KaVo INTRA LUX 68LH/L68/68LD head	15
Sirona T1 Line C200L (Z13 Ø8.7)	16
Midwest E Plus	16
B.A.International BA200LTS	16
Sirona T1 Line C200L (Z13 Ø 9.7)	17
Sirona T1 Line C40 (Ø9.7)	18
Bien Air 1.5(L)	19
NSK Ti-Max Z95(L)	20
NSK Ti-Max X95(L)	21
NSK S-Max M95(L)	22
WH Synea WA-99 LT/A head ø10.0 mm / 3-port spray (old)	23
WH Synea WA-99 LT/A head ø 9.6mm / 5-port spray (new)	23
WH Synea Fusion WG-99 A/LT	24
WH Vision WK-99 LT/LT S	24




KaVo MASTERmatic Lux M25L




KaVo MASTERmatic Lux M05L




KaVo GENTLEpower 25LP/25LPA/25LPR




KaVo EXPERTmatic E25L/E25C



KaVo INTRAmatic 25LH/25LHA/25LHN/25LHC



KaVo INTRA LUX 68LH/L68/68LD head



Sirona T1 Line C200L (Z13 Ø 8.7) Midwest E Plus, B.A.International BA200LTS



Sirona T1 Line C200L (Z13 Ø 9.7)




Sirona T1 Line C40 / C40L, T2 Line A40L / E40



Bien Air 1.5(L)




NSK Ti-Max Z95(L)



NSK Ti-Max X95(L)

i Information



NSK S-Max M95(L)



WH Synea WA-99 LT/A head ø10.0 mm / 3-port spray (old) WH Synea WA-99 LT/A head ø 9.6mm / 5-port spray (new)



WH Synea Fusion WG-99 A/LT WH Vision WK-99 LT/LT S



Parts for High-Speed Handpieces.

On the following pages you will find for different manufacturers and their turbines all suitable spare parts for a complete repair.


KaVo SUPERtorque 625/630/640	28	NSK Ti-Max X500 (NL9000M in the US)
KaVo SUPERtorque 647/649	29	NSK Ti-Max X600 (NL9000 in the US)
KaVo SUPERtorque 650	30	NSK Ti-Max X700 (NL9000T in the US)
KaVo SUPERtorque 655/660	31	Bien Air Bora
KaVo SMARTtorque S609/S619	32	Bien Air Black Pearl
KaVo SMARTtorque S609/S619	32	Bien Air Tornado
KaVo SMARTtorque S609/S619	32	Bien Air Prestige S30-S32
KaVo SMARTtorque Mini S605/S615	33	
KaVo SMARTtorque Mini S605/S615	33	
KaVo SMARTtorque Mini S605/S615	33	
KaVo GENTLEforce 7000	34	
KaVo GENTLEforce 6000	35	
KaVo GENTLEsilence 8000/6500 BR	36	
KaVo GENTLEsilence 6500	37	
KaVo GENTLEmini 4500/5000	38	
KaVo MASTERtorque M9000/M8900	39	
KaVo MASTERtorque Mini M8700	40	
KaVo BELLAtorque 637/Mira Lux 635	41	
KaVo BELLAtorque 639/642/645	42	
KaVo EXPERTtorque E679 L/E680 L	43	
KaVo EXPERTtorque Mini E675 L/E677 L	44	
KaVo COMPACTtorque 636	45	
KaVo POWERtorque 646	45	
Star 430/Solara	46	
Midwest Tradition TB/TL/Pro TBL/Pro TLF	47	
Midwest Tradition Push Button	48	
Midwest Tradition Lever	49	
Midwest Stylus Lite 360/361	50	
Midwest Quiet Air Lever	51	
WH Synea TA 98 (L)	52	
WH Synea TA 98 LED	53	
WH Synea TA 97 LED	54	
WH Fusion TG 98/Synea TA 98 CLED/		
Vision TK-98	55	
WH Vision TK-100	56	
WH Alegra TE 98(L)	57	
Sirona T2/T3 Racer	58	
Sirona T1 Control TC3	59	
Sirona T1 Control CS (SN > 600 000, < 800 000)	60	
Sirona T2 Control / T1 Control (SN > 502.000)	61	
Sirona T2/T3 Boost (SN > 600 000, < 700 000)	62	
Sirona T2/T3 Boost CBM		
(SN > 700 000, < 800 000)	63	
Sirona SiroBoost	64	
NSK Ti-Max Z800	65	
NSK Ti-Max Z900	66	
NSK S-Max M600	67	



KaVo SUPERtorque 625/630/640




KaVo SUPERtorque 647/649



KaVo SUPERtorque 650




KaVo SUPERtorque 655/660




KaVo SMARTtorque S609/S619




KaVo SMARTtorque Mini S605/S615




KaVo GENTLEforce 7000




KaVo GENTLEforce 6000




KaVo GENTLEsilence 8000/6500 BR




KaVo GENTLEsilence 6500




KaVo GENTLEmini 4500/5000



KaVo MASTERtorque M9000/M8900



KaVo MASTERtorque Mini M8700




KaVo BELLAtorque 637/Mira Lux 635




KaVo BELLAtorque 639/642/645



KaVo EXPERTtorque E679 L/E680 L




KaVo EXPERTtorque Mini E675 L/E677 L



KaVo COMPACTtorque 636 **KaVo** POWERtorque 646




Star 430/Solara




	Spare F	Parts	
	Comple	ete rotor/Turbine	1
105 A10549.0	26 39 68	Rotor/Turbine, Lube Free Myonic months g	uarantee
110 A10222.0		Combo 12 months g	uarantee
85 A10719.0		Angular-contact ceramic bearing with flange, Lube Fre	e 🔘
A10492.0		Radial ceramic bearing with flange	
	Back ca	np	2
99 A11026.0		Back cap	
000	O-rings	s/Washers/Misc.	3
88 A10272.0		10-pack set of o-rings	3.1
92 A10330.0		Slinger	3.2

Midwest Tradition TB/TL/Pro TBL/Pro TLF




	Spare F	Parts	
	Comple	ete rotor/Turbine	1
47 A10947.0	66 29 989	premium Rotor/Turbine	months guarantee
102 A10947.1	66 29 989	Rotor/Turbine	months guarantee
109 A10861.0		Combo	months guarantee
84 A10282.0		Radial ceramic bearing	

Midwest Tradition Push Button




	Spare	Parts		
	Compl	ete rotor/Turbine		1
102 A10011.0	790118	Rotor/Turbine	months g	uarantee
109 A10224.0		Combo	nonths g	uarantee
85 A10661.0		Angular-contact ceramic bearing with flang	e	
000	O-ring	s/Washers/Misc.		2
88 A10436.0		10-pack of quad-rings		2.1
92 A10328.0		Slinger		2.2

Midwest Tradition Lever




			Spare P	Parts		
:			Comple	ete rotor/Turbine		1
	102	A10010.0	780118	Rotor/Turbine 12	months g	uarantee
	108	A10370.0		Combo 12	months g	uarantee
	85	A10661.0		Angular-contact ceramic bearing with flange		bur side
	84	A10287.2		Radial ceramic bearing		cancida
	84	A10390.0		Angular-contact ceramic bearing		cap side
		000	O-rings	s/Washers/Misc.		2
	88	A10436.0		10-pack set of quad-rings (blue)		2.1
	88	A10432.0		10-pack of o-rings		2.2
	92	A10328.0		Slinger		2.3

Midwest Stylus Lite 360/361





Midwest Quiet Air Lever




	Spare P	arts		
	Comple	te rotor/Turbine		1
102 A10606.0	484300	Rotor / Turbine	months g	uarantee
108 A10488.0		Combo	months g	uarantee
85 A10661.0		Angular-contact ceramic bearing with flange		
84 A10287.2		Radial ceramic bearing		
84 A10390.0		Angular-contact ceramic bearing		
000	O-rings	/Washers/Misc.		2
88 A10436.0		10-pack of quad-rings (blue)		2.1 (bur side)
88 A10432.0		10-pack of o-rings		2.2 (cap side)
91 A10274.0		10-pack set of wave washers		2.3


WH Synea TA 98 (L)



WH Synea TA 98 LED



WH Synea TA 97 LED




	Spare P	arts		
	Comple	ete rotor/Turbine		1
105 A10458.0	05 30 7200	premium Rotor/Turbine	12 months g	uarantee
105 A10458.1	05 30 7200	Rotor/Turbine	months g	uarantee
111 A10457.0		Combo	months g	uarantee
84 A10288.2	27 70 712	Radial ceramic bearing		
85 A10465.1	27 70 712	Angular-contact ceramic bearing		
000	O-rings	/Washers/Misc.		2
88 A10273.0		10-pack of o-rings		2.1
91 A10274.0		10-pack set of wave washers		2.2
	Light g	uide		3
93 A10715.0	05 28 8700	Glass rod		
	Tools			
94 A10726.0		Back cap key		

WH Fusion TG 98/Synea TA 98 CLED/ Vision TK-98



WH Vision TK-100




	Spare P	Parts		
	Comple	ete rotor/Turbine		1
105 A10806.0	72 72 900	premium Rotor/Turbine	12 months g	uarantee
105 A10806.1	72 72 900	Rotor/Turbine	12 months g	uarantee
111 A10865.0		Combo	12 months g	uarantee
84 A10288.2	27 70 712	Radial ceramic bearing		
85 A10465.1	27 70 712	Angular-contact ceramic bearing		
000	O-rings	s/Washers/Misc.		2
88 A10432.0		10-pack of o-rings		2.1
91 A10274.0		10-pack set of wave washers		2.2

WH Alegra TE 98(L)




Sirona T2/T3 Racer



	Spare P	arts		
	Comple	ete rotor/Turbine		1
104 A10353.0	58 68 307	premium Rotor/Turbine	months g	uarantee
104 A10353.1	58 68 307	Rotor/Turbine	months g	uarantee
109 A10351.0		Combo	months g	uarantee
84 A10288.2	58 68 307	Radial ceramic bearing		
85 A10465.1	58 68 307	Angular-contact ceramic bearing		
000	O-rings	:/Washers/Misc.		2
88 A10432.0		10-pack of o-rings		
	Tools			
94 A10724.0		Back cap key		


Sirona T1 Control TC3



Sirona T1 Control CS (SN > 600 000, < 800 000)



Sirona T2 Control / T1 Control (SN > 502.000)




	Spare P	arts		
	Comple	ete rotor/Turbine		1
104 A10493.0	58 68 299	premium Rotor/Turbine	12 months g	uarantee
104 A10493.1	58 68 299	Rotor/Turbine	12 months g	uarantee
109 A10452.0		Combo	12 months g	uarantee
84 A10288.2	59 94 814	Radial ceramic bearing		
85 A10465.1	59 94 814	Angular-contact ceramic bearing		
000	O-rings	/Washers/Misc.		2
88 A10729.0		10-pack set of o-rings		2.1
91 A10274.0		10-pack set of wave washers		2.2


Sirona T2/T3 Boost (SN>600 000, <700 000)



Sirona T2/T3 Boost CBM $(SN > 700\ 000, < 800\ 000)$



Sirona SiroBoost



	Spare F	Parts		
	Comple	ete rotor/Turbine		1
104 A10843.1	63 34 598	Rotor/Turbine	12 months g	uarantee
109 A10868.0		Combo	•	
84 A10288.2		Radial ceramic bearing		
85 A10465.1		Angular-contact ceramic bearings		
000	O-rings	s/Washers/Misc.		2
89 A10840.0		10-pack set of o-rings		

NSK Ti-Max Z800



NSK Ti-Max Z900



NSK S-Max M600



NSK Ti-Max X500 (NL9000M in the US)




NSK Ti-Max X600 (NL9000 in the US)




NSK Ti-Max X700 (NL9000T in the US)



Bien Air Bora




Bien Air Black Pearl



	Spare P	arts		
	Comple	ete rotor/Turbine		1
106 A10847.0	1500164	Rotor/Turbine	12 months gu	uarantee
111 A10866.0		Combo	12 months gu	uarantee
84 A10358.0		Radial ceramic bearing		
000	O-rings	s/Washers/Misc.		2
88 A10437.0		10-pack set of o-rings		

Bien Air Tornado



	Spare Parts			
	Comple	ete rotor/Turbine		1
106 A10849.0	15 01 968	premium Rotor/Turbine	12 months g	uarantee
106 A10849.1	15 01 968	Rotor/Turbine	12 months g	uarantee
111 A10867.0		Combo	12 months g	uarantee
84 A10288.2		Radial ceramic bearing		
85 A10465.1		Angular-contact ceramic bearing		
000	O-rings	:/Washers/Misc.		2
89 A10840.0		10-pack of o-rings		2.1
91 A10274.0		10-pack set of wave washers		2.2

Bien Air Prestige S30-S32



High-Speed Handpieces

All Products.

On the following pages you will find for different manufacturers and their turbines all suitable spare parts for a complete repair.

Parts/rotors for contra-angle handpieces	78
Bearings	84
O-rings	88
Gaskets	89
LED bulbs	90
Washers	91
Slingers	92
Glass Rod Light Guides	93
Tools/Misc. parts	94
Back caps	98
Turbine couplings	99
Turbines/Rotors for high-speed handpieces	100



Complete heads and head gear cartridges







A10014.0 KaVo 25LP(A)/200XD



A10015.0 KaVo 25LH(A)/25LN



A10016.0 Sirona A/C200L (9.7D), B.A. International 250 LT



A10017.0 Bien Air 1.5L



A10018.0 KaVo E25L/E25C



A10019.0 NSK X95L/M95L



A10090.0 Sirona C200L (8.7D)



A10090.0 B.A. International BA250(L)T, BA200LTS



A10090.0 Sirona A/C200L (8.7D), Midwest E-Plus, B.A. International 200LTS



A10091.0 KaVo M25L



A10092.0 NSK Ti-Max Z95(L)



A10093.0 KaVo M05L



A10852.1 WH WA-99 A



A10853.0 WH WA-99 LT



A10882.0 NSK Ti-Max Z95(L) Cartridge kit



A10902.0 WH Synea Fusion WG-99, Vision WK-99







A10463.0 KaVo 25LP(A) Drive kit



A10519.0 KaVo 25LH(A)/25LN Drive kit



A10525.0 KaVo 25LH(A)/25LN Complete drive set



made in Germany

12 months guarantee

A10550.0 KaVo 25LP(A) Complete drive set



A10564.0 NSK X95L Drive kit



A10565.0 KaVo M25L Drive kit



A10566.0 KaVo M25L Complete drive set



A10574.0 KaVo M05L Drive kit



A10575.0 KaVo M05L Drive kit



A10576.0 KaVo E25L Drive kit



A10625.0 Bien Air 1.5L Drive kit



A10771.0 A10771.0, X95L Head kit



A10774.0 NSK X95L Complete drive set



A10886.0 WH WA-99 LT/A (new) Drive kit



A10887.0 WH WA-99 LT/A (old) Drive kit



A10929.0 Sirona A/C200L (8.7D), Midwest E-Plus, B.A. Intl 200LTS Drive kit









i Information



Gear packs for 1:5 contra-angle intermediate shafts

made in Germany 12 months guarantee



A10620.0

KaVo 25LP(A)/200XD/M25L/ M05L Head side



A10673.0

Bien Air 1.5(L) Head side



A10675.0

KaVo 25LH(A) Head side



A10676.0

KaVo 25LP(A)/M25L/M05L Motor side



A10677.0

KaVo 25LH(A) Motor side



A10770.0

NSK X95L/M95L Head side



A10823.0

Kavo E25 Head side



A10907.0

NSK Z95L Head side



A10920.0

Sirona A/C 200(L) Head side

i Information

two-letter property codes below. Please refer to page 90 for general technical information abou bearings.



Bearings for high-speed / turbine handpieces





A10282.0

CE, RC, TC, AS, ST ID: 3,175 OD: 5.980/6.350 W: 2.78 Germany





A10285.0

CE, RC, TC, AS, XI ID: 3,175 OD: 6,350 W: 2.38/2.78 Germany





myonic myonic

A10287.0

CE, RC, TC, SI ID: 3,175 OD: 6,350 W: 2.38 Germany



A10287.2

CE, RC, TC, AS ID: 3,175 OD: 6,350 W: 2,38 Japan



A10288.0

CE, RC, TC, SI ID: 3,175 OD: 6,350 W: 2,78 Germany



A10288.2

CE, RC, TC, AS
ID: 3,175
OD: 6,350
W: 2,78
Japan



A10291.0

CE, AC, TC, AS, ST, XI ID: 3,175 OD: 5.980/6.350 W: 2.38/2.78 Japan



A10292.0

CE, AC, TC, AS, ST ID: 3,175 OD: 5.980/6.450 W: 2,78 Germany



A10295.0

CE, AC, TC, IS
ID: 3,175
OD: 6,350
W: 2,38
Japan



A10358.0

CE, AC, TC, NS
ID: 3,175
OD: 7,938
W: 3,57
Germany



A10390.0

CE, AC, TC, NS ID: 3,175 OD: 6,350 W: 2,38 Japan

i Information



Bearings for high-speed / turbine handpieces



A10465.1

CE, AC, TC, IS **ID:** 3,175 **OD:** 6,350 **W:** 2,78 **Japan**



A10485.0

CE, AC, TC, IS, ST **ID:** 3,175 **OD:** 5.980/6.350 **W:** 2,78 Japan



A10570.0

CE, RC, TC, AS, FL **ID:** 3,175 **OD:** 7,500 **W:** 2,77 Germany



A10661.0

CE, AC, TC, IS, FL **ID:** 3,175 **OD:** 6,350 **W:** 2,78 **Japan**



A10719.0

CE, AC, PH, IS, FL **ID:** 3,175 **OD:** 6,350 **W:** 2,78 Germany



A10734.0

CE, RC, TC, AS, FL **ID:** 3,175 **OD:** 7,500 **W:** 2,77 Germany



A10830.0

CE, AC, TC, AS, ST, XI **ID:** 3,180 **OD:** 7,08/7,50 **W:** 2,88 Germany



myonic 🍞

A10837.0

CE, AC, TC, IS, FL **ID:** 3,175 **OD:** 6,350/7,900 **W:** 1,95/2,63 Germany



A10838.0

CE, AC, TC, IS **ID:** 3,175 **OD:** 6.350 **W:** 1,95/1,75 Germany



A10839.0

CE, RC, PH, SI, FL **ID:** 3,175 **OD:** 6,350 **W:** 2,38 Germany



A10924.0

CE, RC, TC, IS, ST **ID:** 3,175 **OD:** 6,350 **W:** 2,35/2,45 Germany



i Information



Bearings for contra-angle / attachment handpieces



A10325.0 CE, RC, TC, AS **ID:** 4,000 **OD:** 7,000 **W:** 2,00

Germany



A10326.0 CE, RC, TC, AS, XI **ID:** 3,175 **OD:** 6,350 **W:** 1.90/2.70 Germany



A10335.0 CE, RC, TC, NS **ID:** 3,175 **OD:** 7,000 **W:** 2,00 Germany



A10336.0 CE, RC, TC, NS **ID:** 3,175 **OD:** 6,000 **W:** 2,00 Germany



A10371.0 SS, RC, TC, AS, XI **ID:** 3,000 **OD:** 7,000 **W:** 2.35/3.35 Germany



A10372.0 SS, RC, TC, AS **ID:** 4,000 **OD:** 7,000 **W:** 2,00 Germany



A10396.0 SS, RC, TC, AS **ID:** 5,000 **OD:** 8,000 **W:** 3,00 Germany



A10464.0 SS, RC, TC, NS **ID:** 2,000 **OD:** 5,000 **W:** 1,50 Germany



A10491.0 SS, RC, TC, NS **ID:** 3,000 **OD:** 7,000 **W:** 3,00 Germany



A10621.0 CE, RC, TC, NS **ID:** 4,000 **OD:** 7,000 **W:** 2,00 Germany



A10663.0 SS, RC, TC, AS **ID:** 2,350 **OD:** 5,000 **W:** 1,50 Germany



A10666.0 CE, RC, TC, AS **ID:** 4,000 **OD:** 7,000 **W:** 1,60/2,00 Germany

i Information



Bearings for contra-angle / attachment handpieces



A10667.0 CE, RC, TC, NS, XI **ID:** 3,175 **OD:** 6,350 **W:** 1.92/2.78 Germany



A10765.0 SS, RC, TC, NS **ID:** 2,380 **OD:** 4,763 **W:** 1,59 Germany



A10781.0 SS, RC, TC, NS **ID:** 2,000 **OD:** 4,700 **W:** 6,0 Germany



A10821.0 CE, RC, TC, NS, XI **ID:** 3,175 **OD:** 5,800 **W:** 1,42/2,28 Germany



A10822.0 CE, RC, TC, AS **ID:** 3,500 **OD:** 6,500 **W:** 2,00 Germany



A10834.0 CE, AC, TC, SI **ID:** 3,800 **OD:** 6,600 **W:** 1,75 Thailand



A10925.0 CE, RC, TC, NS **ID:** 2,358 **OD:** 5,500 **W:** 2,00 Germany



A10926.0 CE, RC, TC, NS **ID:** 3,400 **OD:** 6,500 **W:** 1,60 Germany

O-rings



A10271.0

10-pack set of KaVo o-rings



A10272.0

10-pack set of o-rings for Star Solara turbines



A10273.0

10-pack set of o-rings for WH, Sirona T2/T3 Racer



A10432.0

10-pack set of o-rings 6.0 x 0.85



A10436.0

10-pack blue Quadrings Midwest



A10437.0

10-pack set of o-rings for Bien Air Bora



A10653.0

10-pack set of o-rings for head rotor Bien Air 1.5L



A10669.0

O-ring set for MULTIflexcompatible coupling



A10672.0

10-pack set of o-rings for NSK X600(L)/NL9000/MW Stylus



A10704.0

10-pack set of o-rings 6.2 x 0.80



A10729.0

10-pack set of o-rings for Sirona T2 Control (6.3 x 0.8 mm)



A10745.0

10-pack O-ring d 4.9 x 0.3 mm for X95L/M95L/Z95L



A10746.0

10-pack O-ring d 6.9 x 0.3 mm for canister X95L/M95L



A10763.0

10-pack O-ring for intermediate shaft X95L/M95L



A10797.0

10-pack O-ring for spray tube, rear (motor) side KaVo 25LP/E25L



A10798.0

10-pack O-ring for spray tube, front (head) side KaVo M25L/25LP/E25L

O-rings/Gaskets



O-rings



A10799.0

10-pack O-ring for spray tube, rear (motor) side KaVo M25L



A10800.0

2-pack O-Ring for idler gear KaVo 25LP(A)/ M25L and intermediate shaft 25LH(A)



A10840.0

10-pack set of o-rings



A10919.0

O-ring set for NSK coupling



A10921.0

10-pack set of o-rings for NSK M600 (5.25x0.3mm)



A10922.0

10-pack set of o-rings for NSK M600 (6.90x0.3mm)



Gaskets



A10766.0

2-pack of gaskets for NSK Ti/S-Max X95(L)/M95(L)



A10872.0

Gasket for ISO 6-hole connection, turbine/coupling

LED bulbs

coupling + motor



Washers



A10274.0

10-pack set of wave washer for KaVo and WH turbines



A10327.0

10-pack set for KaVo turbines 6000/6500/4500/5000/E675/ E677/E679/680



A10347.0

10-pack set of wave washers for cartridge KaVo 25LP(A)/E25L, NSK Z95L



A10733.0

10-pack set of wave washers and distance washer for Bien Air 1.5L



A10824.0

10-pack wave washer for KaVo M25L



A10828.0

10-pack set of WH Alegra TE 98(L) wave washers



A10875.0

10-pack set of Sirona T2/T3 Boost wave washers



A10893.0

10-pack set of wave washers for cartridge NSK X95L/M95L



A10923.0

10-pack set set of wave washers for Sirona/Midwest special stepped bearing

Slingers



Slingers



A10364.0 Slinger for KaVo 635/637



A10316.0 Slinger for KaVo 650



A10328.0 Slinger for Midwest Tradition Lever/Push-Button



A10329.0 Slinger for KaVo 647/649



A10330.0 Slinger for Star 430/Solara



A10364.0 Slinger for KaVo 635/637



A10532.0 Slinger for KaVo 8000/6500, NSK X600



A10581.0 Slinger for NSK Ti-Max Z800



A10688.0 Slinger for Midwest Stylus



A10652.0 Distance washer head rotor Bien Air 1.5L



A10894.0 5-pack set of slinger washer for cartridge NSK X95L/M95L



A10895.0 5-pack set of slinger bushing for cartridge NSK X95L/M95L/Z95L



A10896.0 5-pack set of slinger washer for cartridge NSK Z95L




Slinger set for cartridge NSK Z95L

Glass Rod Light Guides

A10739.0

Glass rod for KaVo 7LH/20LH/29LH/2068LH A10740.0

Glass rod for KaVo 20LP / E20L



A10773.0

Glass rod for WH WK-93/99 LT

A10741.0

Glass rod for KaVo 24LN

Glass Rod Light Guides/Tools/Misc. parts

For high speed turbines and attachments



A10784.0 Glass rod for NSK Z-Series contra-angles



A10785.0 Glass rod for NSK M-Series contra-angles



A10820.0 Glass rod for KaVo SONICflex 2003 / 2008L

Back caps key



A10654.0 Back cap Key KaVo 655/660/6000/6500



A10655.0 Back cap key KaVo 630/640/642



A10720.0 Back cap key for Bien Air Bora



A10721.0 Back cap key for KaVo 5000 / M9000



A10723.0 Back cap key for Sirona T1 Control / T1 Mini



A10724.0 Back cap key for Sirona T2/T3 Racer



A10725.0 Back cap key for WH Synea TA 98 / TA 96



A10726.0 Back cap key for WH Synea TA 97 LED/TA 98 LED

Tools/Misc. parts

More information on the technical specifications of our oil, including the adapters needed for different handpiece types can be found at

https://www.adensys.net/product-category

i Information



Back caps key



A10727.0Back cap key for
NSK Ti-Max X500/X600/X700



A10728.0Back cap key for KaVo 635/637



A10730.0 Key for crown cap on KaVo 25LP(A)



A10726.0 Back cap key for WH TE-97/TE-98



A10736.0 Back cap key for NSK X95L/M95L



A10778.0 Back cap key for NSK M600/Pico



A10779.0Back cap key for WH TK 94/TK 100



A10888.0 Back cap key for NSK Z900



A10909.0 Key for crown cap on KaVo M25L



A10910.0Back cap key for Sirona T1, T2
Control (600 000 < SN < 800 000)



A10913.0Back cap key for Sirona T1-T3
Boost (600 000 < SN < 800 000)



A10927.0 Back cap key and reamer for NSK Z800

Tools/Miscellaneous parts

i Information



Other tools / parts / maintenance equipment



A10299.0 Retention force checker (for 1.6 mm bur)



A10668.0 Digital one-channel manometer for drive air



A10678.0 Removal tool for front gear of interm. shaft, D4.0 mm (ex: X95L, Z95L, Bien Air 1.5L)



A10731.0 Ejector pin for push button KaVo 25LP(A)



A10744.0 ADENSYS precision tweezers 105mm stainless steel, non-magnetic



A10783.0 Water filter for KaVo 6000/6500/7000/8000



A10790.0 ADENSYS handpiece maintenance spray 500ml



A10793.0 Spray adapter for MULTIflex/KaVo



A10794.0 Spray adapter for Sirona



A10795.0 Spray adapter for NSK



A10796.0 Spray adapter for WH



A10809.0 ADENSYS handpiece maintenance spray 250ml



A10813.0 6-pack ADENSYS handpiece maintenance spray 500ml



A10878.0 3-in-1 spray adapter for contra-angles



A10891.0 Jig to hold intermediate shaft during assembly / dis-assembly



A10898.0 Tool for o-ring assembly inside canister NSK

Tools/Miscellaneous parts

Information

More information on the technical specific

handpiece types can be found at https://www.adensys.net/product-category/



Other tools / parts / maintenance equipment



A10899.0

Tool to remove rear bearing of NSK X95L/M95L



A10900.0

Tool to push actuator side of turbine spindles



A10901.0

Fixture to remove front bearing + slingers, NSK X95L/M95L/Z95L



A10903.0

Fixture to help remove head housing KaVo 25LP(A)



A10904.0

Fixture to help remove head housing KaVo E25L



A10905.0

Fixture to help remove head housing KaVo M25L/M05L



A10928.0

Jig to remove canister NSK X95L/ M95L

Back caps



Back caps for high-speed / turbine handpieces



A10035.0 Back cap for NSK Ti-Max X500



A10742.0 Complete back cap for KaVo E25L



A10914.0 Complete back cap for NSK Ti-Max X700



A10915.0 Complete back cap for Bien Air Black Pearl



A10917.0 Complete back cap for NSK Ti-Max X95(L)



A11013.0 Complete back cap for KaVo 640



A11014.0 Complete back cap for Bien Air Bora



A11015.0 Complete back cap for KaVo 635 / 637



A11016.0 Complete back cap for KaVo 647 / 649 /650



A11017.0 Complete back cap for KaVo 6500 / 8000 / M9000



A11018.0 Complete back cap for KaVo 655 / 660 / 6000 / 7000



A11019.0 Complete back cap for WH TA-98



A11020.0 Complete back cap for WH TA-98 LED



A11021.0 Complete back cap for . KaVo E679 / E680



A11022.0 Complete back cap for KaVo 639 / 642 / 645



A11023.0 Complete back cap for WH TE-98 RM

Back caps/Turbine couplings



Back caps for high-speed / turbine handpieces



A11024.0 Complete back cap for NSK X600L



A11025.0 Complete back cap for NSK M600L



A11026.0 Complete back cap for Star 430



A11027.0 Complete back cap for KaVo 4500 / 5000 / E675 / E677



A11029.0 Crown cap for KaVo 25LP/25LPA/25LPR



A11030.0 Complete back cap assembly for KaVo M25L



A11034.0 Complete back cap for NSK Ti-Max Z900



A11035.0 Complete back cap for NSK Ti-Max Z800

Turbine couplings



A10660.0 MULITflex (KaVo)-style 6-pin LED coupling



A10679.0 NSK compatible 6-pin LED coupling



A10780.0 Sirona 6-pin LED coupling for turbine / high-speed handpiece



A10836.0 WH RQ-24 compatible 6-pin coupling



Complete rotors/turbines









A10007.0KaVo 6000 premium



A10007.1 KaVo 6000



A10008.0KaVo 6500B premium



A10008.1 KaVo 6500B



A10250.0 KaVo 630/640 premium



A10250.1 KaVo 630/640



A10251.0 KaVo 650 premium



A10251.1 KaVo 650



A10252.0 KaVo 655/660 premium



A10252.1 KaVo 655/660



A10253.0 KaVo 7000 premium



A10253.1 KaVo 7000



A10254.0KaVo 8000/6500 BR premium



A10254.1 KaVo 8000/6500 BR



A10255.0 KaVo 635/637 premium



A10255.1 KaVo 635



Complete rotors/turbines







A10259.0 KaVo 642/645 premium



A10259.1 KaVo 642/645



A10300.0 KaVo 647/649 premium



A10300.1 KaVo 647/649



A10345.0 KaVo 4500/5000 premium



A10349.0 KaVo 635/637



A10349.1 KaVo 635/637



A10533.0 KaVo 636/646 premium



A10533.1 KaVo 636/646



A10535.0 KaVo E679/E680 premium



A10535.1 KaVo E679/E680



A10540.0 KaVo M8900/M9000 premium



A10540.1 KaVo M8900/M9000



A10560.0 KaVo A675/677 premium



A10560.1 KaVo E675/E677



A10845.1 KaVo S609/S619



Complete rotors/turbines









A10844.1

KaVo Mini S605/S615





A10846.0

KaVo Mini M8700 premium



A10846.1

KaVo Mini M8700

Midwest



A10010.0

Midwest Tradition Lever



A10011.0

Midwest Tradition Push-Button



A10431.0

Midwest Stylus 360/361



A10606.0

Midwest Quiet Air Lever



A10947.0

MidwestTraditionTB/TL/ProTBL/ ProTLF premium (2017-2019)



A10947.1

Midwest Tradition TB/TL/Pro TBL/Pro TLF (2017-2019)



Complete rotors/turbines









A10394.0 NSK TiMax X600 (NL9000)



A10556.0 NSK S-Max M600



A10571.0 NSK Ti-Max Z900 (original-style bearings)



A10571.1 NSK Ti-Max Z900 (alternative-style bearings)



A10580.0 NSK Ti-Max Z800



A10801.0 NSK TiMax X700 (NL9000T) premium



A10801.1 NSK TiMax X700 (NL9000T)



A10850.0 NSK Ti-Max X500



Complete rotors/turbines





Sirona



A10353.0 Sirona T2/T3 Racer premium



A10353.1 Sirona T2/T3 Racer



A10456.0 Sirona T1 Control (TC3) premium



A10456.1 Sirona T1 Control (TC3)



A10493.0 Sirona T2 Control premium



A10493.1 Sirona T2 Control



A10563.0 Sirona T1 Control CS



A10805.1Sirona T2/T3 Boost (SN > 600.000)



A10854.1 Sirona T2/T3 Boost CBM (SN > 700.000)



A10843.1 Sirona SiroBoost



Complete rotors/turbines









A10337.0 WH TA-98 premium



A10337.1 WH TA-98



A10558.0 WH Alegra TE-98(L)



A10359.0 WH TA-98 LED premium



A10359.1 WH TA-98 LED



A10360.0 WH TA-96 premium



A10360.1 WH TA-96



A10601.0 WH TG-98/TA-98/TK-98L CLED premium



A10601.1WH TG-98/TA-98/TK-98L
CLED



A10458.0 WH TA-97 LED premium



A10458.1 WH TA-97 LED



A10806.0 WH Synea Vision TK-100 premium



A10806.1 WH Synea Vision TK-100

made in Germany

12 months guarantee










BA 675/680/685 Ultimate









Turbines/Rotors for high-speed handpieces



Turbines/Rotors for high-speed handpieces



Turbines/Rotors for high-speed handpieces











Bien Air Bora





Dental bearing technology.

For all rotating dental equipment, bearings are a critical determinant of the performance and overall quality. Given that most dental equipment undergoes reprocessing between patients (including disinfecting and/or sterilizing), bearings undergo a lot of thermal and chemical stresses. For this reason, it is extremely important that high quality bearings are used in order to ensure optimal product reliability.

At ADENSYS, we work with bearings from a number of different premium-quality manufacturers based in either Japan or Germany. One thing is clear – bearings are all different! One can never just compare, for example, $4 \times 7 \times 2$ mm bearings without context. It is important to understand the different parameters and mechanical properties when installing and using bearings to service dental equipment.

We know how confusing it can be to understand the different factors that separate bearings from one another. Please consult the following guide to help you through the dizzying world of dental bearing technology..

Steel Versus Ceramic

If you are involved in the repair of dental handpieces, you will surely have been confronted by the question of "steel versus ceramic" bearings. You may have even heard that ceramic are better. But what does this mean? First off, we need to be clear about nomenclature: when we talk about ceramic bearings, these are in fact hybrid bearings. They are built with steel inner and outer races, but it is the balls that are ceramic. The ceramic balls are typically silicon nitride, or an equivalent material.

Ceramic bearings do, indeed, offer two key benefits over steel bearings:

- 1. Ceramic balls are rounder than steel balls. They also have a smoother surface and more uniform size. These two factors lead to lower overall rolling resistance, which in turn lead to lower turning friction for the bearing. Lower friction means that less torque is needed to turn the bearing, and more input energy is transferred to the spindle.
- 2. Ceramic balls are harder than steel balls, and the ceramic material will never corrode. This improves durability and will last to longer lifetime overall.

These two factors mean that, if all other factors are equal, a ceramic bearing will outlast and outperform a steel bearing at all times. For this reason, ADENSYS carry only ceramic bearings in their product portfolio for all high speed applications (i.e. >150K rpm).

ABEC or ISO quality levels

Bearings are classified in quality levels, which correspond to the precision of the individual parts that are being used to build the final product. In the US, these quality levels are classified according to the ABEC scale. In Europe, this quality is classified on the ISO P-scale. Looking at the ABEC scale, the higher the ABEC classification, the more precise the individual components. For the P-scale, it is the opposite; the lower the P-classification, the higher the part precision. Generally, for dental handpiece applications, we suggest the use of ABEC 5 / P5 or better. Although better (i.e. more precise) classifications are generally better for bearing lifetime and performance, the prices also increase with every step up. We suggest that it does not make sense to install any bearings with classifications above ABEC 9 or P4, as the precision improvements above these levels will have negligibe effect on performance and lifetime, despite large cost increases.

At ADENSYS, we generally use ABEC 7/9 or P4P bearing classifications for our high-speed applications (i.e. for turbines or the head gear rotors in 1:5 contra-angles). These are hybrid quality levels (with some components at ABEC 7 / P4 quality levels and some components at even higher levels).

Please be sure that, whenever you compare different bearings, you know the quality classification! You should also check that marketed classifications are legitimate. It is common to see companies selling dental bearings that have been promoted as "ABEC 9", for example, but when you ask for documentation to prove this, none exists.

Cage materials: torlon and phenolic

The bearing cage is a part of the bearing that resides between the inner and outer races. Its main function is to keep the space between balls as the bearing rotates. The interface between the balls and the cage account for a large part of the friction inside a bearing and the cage material is therefore selected to introduce as little friction as possible within the bearing. At the same time, the cage need to be able to withstand the high temperature and corrosive loads that come from the reprocessing of dental handpieces. When bearings fail within a dental handpiece or motor, the cause of the failure can be traced to the failure of the cage in a huge majority of cases. Generally, there are two materials that are found in the majority of high-speed dental bearings:

Phenolic



This cage type is actually a cloth material, which is set in place using a phenolic resin. It can be identified by its brownish colour. A phenolic cage is desirable for use in dental applications for the following reasons:

- · Low coefficient of friction
- Ability to absorb lubricants
- High strength-to-weight ratio
- · Easily machined

The one drawback of the phenolic cage is that it does not deal well with repeated thermal loads (like from a autoclave/sterilizer). For this reason, engineers have moved to increasingly using different cage materials for dental handpieces.



Torlon

This cage type is a polyamide resin that is strengthened with graphite. The graphite addition gives it its black colour. Torlon has many of the same benefits of phenolic, but it has improved temperature resistance and excellent chemical resistance. The one drawback of the torlon cage is that it does not absorb lubrication as well as phenolic, meaning that it is more sensitive to missing lubrication.

It should be noted that bearings cages from both of these materials can be produced through machining or a moulding process. Machined cages will generally offer better stability and overall mechanical strength but machined cages are generally more costly. All of this is to say that it is not just important to know the cage material, but it is also important to know how the bearings was produced!


Radial versus Angular-Contact

Generally, in dental handpiece applications, you will encounter two main bearing geometry/ construction principles: 1) deep-groove ball bearings (DGBB), also known as radial bearings, or 2) angular contact ball bearings (ACBB). Each of these bearing types has distinct advantages and disadvantages. As you will see in this catalogue, in the section on dental bearings, ADENSYS carries both DGBB and ACBB in inventory. This way, we can address customers who express a preference for either of these bearing types. Still, it is important to understand the differences.



DGBB

Radial bearings are produced with symmetrical grooves in the inner and outer races in which the balls travel as the bearing rotates. They can therefore accept axial loads in both directions. Related to dental handpiece repair, this means that there is no risk of the bearing falling apart when it is assembled or disassembled from the handpiece.




However, given the building principle of these bearings, the cage is assembled as a final step and must be fitted in place from one side. For this reason, radial bearings use a so-called snap-in cage. The disadvantage of a snap-in cage is that it has less mechanical stability than a full cage (see ACBB). This makes a radial bearing inherently less robust, in terms of bearing lifetime, than an angular contact bearing.



ACBB

Angular contact bearings have an asymmetrical geometry, given that the grooves in the inner and outer bearing races are machined with a particular contact angle. Radial and axial loads are designed to be transferred from one ring to the other, through balls, along a specific contact line. With this design, there is nonstop contact between the balls and both bearing raceways. This allows for much higher speed ratings. The contact line is determined by a taper, which is machined into the contact surface of the inner or outer race of the bearing.



A massive benefit of the ACBB is that the build principle call for the tapered element to be assembled as last step (depending on the design type, this means either the inner or outer race will be assembled last). Given that the bearing cage is not assembled as a last step, it is possible to use a full bearing cage design. The full bearing design is much more robust thatn the snap-in cage used in radial bearings.

The disadvantage of an ACBB is that the bearing cannot support axial loads in both directions. With regards to dental handpiece repair, this means that, although an ACBB can support huge load in the operational axial direction, the bearing can fall apart when it is pulled in the opposite direction. This commonly happens if a rotor or turbine assembly is pulled out of a head housing prior or during a repair.

As explained above, the ACBB has huge advantages in terms of overall performance and reliability. For this reason, most dental handpiece manufacturers use angular contact bearings in new handpieces. However, the weakness of the ACBB to support axial loads in both directions means that many repair professionals choose to use radial bearings for dental handpiece repair.

Conclusion

As we have written a few times in this text, there are a huge amount of mechanical and construction properties that can differentiate one bearing from another, even if they have the same general outer dimensions. To be clear, this text only highlights some of the most relevant properties for handpiece repair technicians. There are still many more properties that could be analysed (internal clearances, shield construction, lubrication and much more...).

The most important message that we want to bring across is that bearings should not be compared on an apples-to-apples basis until all of these factors have been understood. Next time you are offered a great price for a dental bearing, we hope you will use this guide to understand more about the real value. Not just price!

Please contact us if we can support you in finding the right bearing for your application.

The ADENSYS Quality Standard.

As the two founders of ADENSYS, we have combined experience of **more than 40 years** in the development and production of dental equipment. The formative years of this product development experience were spent at **world-renowned manufacturers** in the German-speaking dental world, and the years in our pre-ADENSYS careers helped craft our **vision** for **value**, **quality and reliability**. We realize that the dental market, like any other industry, is competitive.

Therefore, we know we have to bring **cost-effective solutions** to the market. But, at the same time, we believe that customers are best served when they **receive products and services** that go above and beyond. We are not satisfied to deliver just another "me-too" product at a **better price**.

It is **our goal** to exceed **customer expectations** in terms of the **quality and performance** of our products and services. If we can deliver this at **competitive price**, we've accomplished our mission!

Given this, ADENSYS strives to be the **highest quality provider of parts** and services regarding the **manufacturing and service** of dental handpieces.

The ADENSYS Mission

It is easy to write about the high quality we strive to deliver but it is a whole other thing to deliver on this every day, for every transaction or product. To start with, we live the mantra indicated above in all of our processes and dealings, and we extend this philosophy to all of our stakeholders: our shareholders, our suppliers, our employees and our customers. Additionally, our philosophy has to extend to all of our business processes:


- 1. We use processes to ensure quality in our product and service development
- 2. We build checks and assurances prior to product or service launch
- 3. We monitor results in the markets in order to sense issues that need immediate addressing
- 4. We define and implement corrective actions

This is dedicated implementation of the PDCA cycle, which was developed by the renowned engineer and statistician Walter Shewhart in 1929:



ADENSYS Quality – materials and construction





Turbine Dr.



Turbine Dr. - Get this new smartphone app from ADENSYS!

The Turbine Dr. App allows you to reliably check the optimal performance of your dental turbines for improved productivity and more safety. The app is now available for Android and iPhone for download.

With the Turbine Dr. App you will immediately know the current operating status of your turbines - **quick and easy!**

ADENSYS recommends the routine performance check of air driven turbines. More and more dentists realize that it is both cheaper and more efficient to ensure the optimal performance of dental tools in advance to avoid surprises during use.

Turbine Dr. offers a revolutionary diagnostic technique for rapid performance testing of their turbines, essential for comfortable working and successful treatments. Use the app even after a self repair to ensure ultimate efficiency.

Diagnosis Mode

With Turbine Dr. in diagnostic mode you can check the current operating condition of your air driven turbines. The app uses the sound signal, checks the complex interaction between the drive air pressure and the turbine speed, and detects faults with great accuracy.

RPM Mode

The speed check mode "RPM" is for accurate checking of the efficiency of your turbines and for checking a self repair. Just let your turbine run at full speed next to the smartphone and the Turbine Dr. App will immediately display the measured speed.



Liability

Any liability – for whatever legal reasons – for loss or damage accruing that are not due to the delivered products shall only be accepted by us if such loss or damage is caused by negligence of a substantial contract obligation or by gross negligence, intent or malice by us or if we have guaranteed the absence of a defect. If we shall be liable for the negligence of a substantial contract obligation, our liability shall be limited to such loss or damage which might typically be predicted.

Moreover we shall be liable for claims according to the Product Liability Act and for damages due to culpable injury to life, body or health.

Any other claims are excluded.

Trademarks

All brands and models that are mentioned in ADENSYS marketing or sales material are registered trade names. No implication is made that any of our products are original products or that any of the referenced companies endorse our products. This clause applies to the following names:

KaVo, Star Dental, Midwest, NSK, GENTLEsilence, GENTLEforce, MULTIflex, WH, Sirona, Bien Air, Sirona, MASTERmatic, GENTLEpower, EXPERTmatic, INTRAmatic, INTRA LUX, Ti- Max, S-Max, SUPERtorque, MASTERtorque, BELLAtorque, EXPERTtorque, GENTLEmini, COMPACTtorque, POWERtorque, Solara, Stylus, Quiet Air, Synea, Alegra, Bora

ADENSYS AG

Messinastrasse 30 9495 Triesen Liechtenstein www.adensys.net info@adensys.net +423 392 6566

ADENSYS GmbH

Hauptstrasse 41 87740 Buxheim Germany www.adensys.net info@adensys.net +49 8331 98 46 897